


M A X I C O M F O R T TM S E R I E S




Day Dreamer
PR632
shown in Coffee Bean



Hazelnut



Calypso



Shiraz



Coffee Bean*



Bridle*

* upcharge applies



PR632

The DayDreamer™ is a medium sized power lift recliner offering our exclusive PowerPillow™ articulating headrest and The Cloud™ seating system, featuring our custom bucket seat design. Recline, relax, and adjust your head position with the touch of a button.





	Day Dreamer
Model Number	PR632
Size	Medium
Weight Capacity	375 Lb.
Recline Positions	MaxiComfort
Standard Fabrics	5
Fabric Options	ALL
Chaise Seat	YES
Coil Springs	No
Back Type	Seam
Head Rest Cover	NO
Arm Covers	NO
Battery Back Up	YES
Overall Width	38.5"
Overall Height	45"
Width Between Arms	22.5"
Floor to Top of Seat	20"
Seat Depth	20.5"
Seat to Top of Back	31"
Seat to Top of Arm	6.5"
Seat Angle Raised	15°
Floor to Top of Foot Rest Reclined	24"
Floor to Top of Back Reclined	23.5"
Back Angle Closed	105°
Back Angle Reclined	165°
Foot Rest Angle Closed	90°
Foot Rest Angle Opened	170°
Distance Required from Wall Reclined	22"
Extension reclined from Back of Seat to end of Foot Rest	42.5"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 Yrs.
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 Yrs.
Warranty on Electrical Labor	1 Yr.
Removable Back Option Available	NO
Smart Tek System Standard	YES
Customizable	NO
Left Hand Control Option Available	YES
Extra Pocket Option Available	Standard
Heat & Massage Option Available	YES
Gel Memory Foam Option Available	NO
Power Pillow Option Available	Standard
AutoFlex Option Available	NO
Foot Rest Extension Option Available	NO
Straight Lift Option Available	NO
User Height (suggested)	5'4" - 5'10"

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Power Articulation Head Forward



Power Articulation Head



Power Articulation Head Back



Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with
UL 73 and
CSA C22.2 No. 68

